

## UNISTATION<sup>TM</sup>

2<sup>nd</sup> Version

The simplest & easiest way to various cellular therapies

Unique all-in-one solution







## **UNISTATION<sup>™</sup>**

2<sup>nd</sup> Version

The simplest & easiest way to various cellular therapies

Unique all-in-one solution





# 

#### **Background of Development**

- Zuk method
- CAL

#### **Detailed Information**

- Composition
- Specification
- Consumables & Accessories

#### Protocols

- ADSC isolation
- Cell counting

#### **Comparison & References**

- Comparison of other devices
- References



## Upgrades on UNISTATION™ 2<sup>nd</sup> Version

The simplest & easiest way to various cellular therapies

Unique all-in-one solution



#### Upgrades on UNISTATION 2<sup>nd</sup> version





Plain design

Modern and state of the art design

#### Upgrades on UNISTATION 2<sup>nd</sup> version



#### UniStation 1<sup>st</sup> Version (2013)



Some parts not autoclavable (eg. Tuberack and shaking plate)

#### UNISTATION<sup>™</sup> 2<sup>nd</sup> Version





<u>Autoclavable</u> All parts are autoclavable



#### Button push style keypad



<u>Touch screen keypad</u> No button push/ Modern style



Air heater type incubator Normal incubation method *Silicon heater type incubator* Safer and faster incubation method

#### Copyright © 2016 NeoGenesis Co., Ltd.

#### Upgrades on UNISTATION 2<sup>nd</sup> version



#### UniStation 1<sup>st</sup> Version (2013)



#### **UNISTATION<sup>™</sup> 2<sup>nd</sup> Version**

- U.S. Department of Health & Human Se	rvices					a A A
FDA U.S. Food an Protecting and	<b>nd Drug Admin</b> Promoting <i>Your</i> He	istration ealth		A to Z In Most Pa	dex   FollowFDA   En pular Searches	Español SEARCH
Home Food Drugs Medical Dev	vices Radiation-Emitting	Products Vaccines, E	Blood & Biologic	cs Animal & Vet	erinary Cosmetics	Tobacco Products
Establishment Regis FDA Home O Medical Devices O	tration & Dev	vice Listing				a 🛛 🖂
	New Search			Back To Search R	esults	
	Proprietary Name: Classification Name: Device Class: Regulation Number: Medical Specialty: Registered Establishment Name: Owner/Operator: Owner/Operator Number: Establishment Operations	CENTRICite 1000 and CENTRIFUCES (MICR LINICAL USE JQC 1 882.2050 Clinical Chemistry <u>NEOGENESIS</u> 10048015 : specification Develope	Accessories IO, ULTRA, REF	RIGERATED) FOR	2	
Page Last Updated: 12/01/2014 Note: If you need help accessing information	on in different file formats, s	ee Instructions for Down	FOIA No	s and Players. • Fear Act   Site	e Map   Transparen	cy   Website Policies
U.S. Food and Drug Administra 1993 New Hampshire Avenue Silver Spring, MD 20993 Ph. 1888-NP-CPA (1-888-463-6332) Emai FDA Statoz, Statoz, Stato For Government   For Press	tion Count Advis Regu Safet Emm Inter Inter Safet Safet Inter Inte	sination Products sory Committees one & Research Latery Information Y gency Preparedness a Events is Events in grand Confining Education compliance & Local Officials umors by n Professionals tructive	n	🦧 U.S. Depi	artment of <b>Health &amp; H</b>	uman Services
	FDA (	US) re	egis	tered	<u>d</u>	

FDA, CE, KFDA

FDA (US) not registered CE, KFDA

#### Introduction





## UNISTATION™ is, All-in-one Device

for diverse autologous cell therapies

UNISTATION<sup>™</sup> is a specially designed medical device to satisfy diverse demands for various autologous cell therapies with only one device.

<u>Youtube Link of UNISTATION Introduction</u> http://youtu.be/-Jy3t4y1ny0

<u>Youku Link of UNISTATION Introduction (Chinese)</u> http://v.youku.com/v\_show/id\_XODI0MjUzODk2.html

#### Application areas of various treatments















#### What is Stem Cell?

- Stem cells are multi potent cells and its multi-diversification capability is considered as the core future of medical fields demanding various cell therapies for a number of incurable diseases.
- There are various types of stem cells, embryonic stem cells, umbilical stem cells, and adult stem cells and adult stem cells includes adipose derived stem cells, peripheral blood stem cells, bone marrow stem cells, and etc.
- However, embryonic stem cells have controversial ethical problems and possible harvesting period of umbilical stem cells are very limited to the time of delivery, so the use of them from human sources are extremely prohibited or limited.
- Thus, the need for harvesting a large amount of healthy adult stem cells is increasing.

#### What is **ADSC** (Adipose Derived Stem Cell)?

- Adipose tissue is one of the richest sources of MSCs. When compared to bone marrow, there is more than 500 times more stem cells in 1 gram of fat when compared to 1 gram of aspirated bone marrow. Adipose stem cells are currently actively being researched in clinical trials for treatment in a variety of diseases.
- In 2001, Dr. Patricia Zuk first published a thesis explaining how stem cells from adipose tissue were successfully obtained and cultured by their research team. (Tissue Engineering. 7:212-228, 2001)
- Due to the spreads of the publication, adipose tissue started to be known as the best source of stem cells among scholars and doctors due to the easy harvesting technique that Dr Zuk showed in his thesis and due to the obtainability in large amounts as well.
   *Reference: Tissue Engineering. 7:212-228, 2001*



TISSUE ENGINEERING Volume 7, Number 2, 2001 Mary Ann Liebert, Inc.

#### Multilineage Cells from Human Adipose Tissue: Implications for Cell-Based Therapies

PATRICIA A. ZUK, Ph.D.,<sup>1,2</sup> MIN ZHU, M.D.,<sup>1,2</sup> HIROSHI MIZUNO, M.D.,<sup>2</sup> JERRY HUANG, B.S.,<sup>2</sup> J. WILLIAM FUTRELL, M.D.,<sup>3</sup> ADAM J. KATZ, M.D.,<sup>3</sup> PROSPER BENHAIM, M.D.,<sup>2</sup> H. PETER LORENZ, M.D.,<sup>2</sup> and MARC H. HEDRICK, M.D.<sup>2</sup>

#### ABSTRACT

Future cell-based therapies such as tissue engineering will benefit from a source of autolo gous pluripotent stem cells. For mesodermal tissue engineering, one such source of cells is the bone marrow stroma. The bone marrow compartment contains several cell populations, including mesenchymal stem cells (MSCs) that are capable of differentiating into adipogenic, osteogenic, chondrogenic, and myogenic cells. However, autologous bone marrow procurement has potential limitations. An alternate source of autologous adult stem cells that is obtainable in large quantities, under local anesthesia, with minimal discomfort would be advantageous. In this study, we determined if a population of stem cells could be isolated from human adipose tissue. Human adipose tissue, obtained by suction-assisted lipectomy (i.e., liposuction), was processed to obtain a fibroblast-like population of cells or a processed lipoaspirate (PLA). These PLA cells can be maintained in vitro for extended periods with stable population doubling and low levels of senescence. Immunofluorescence and flow cvtometry show that the majority of PLA cells are of mesodermal or mesenchymal origin with low levels of contaminating pericytes, endothelial cells, and smooth muscle cells. Finally, PLA cells differentiate in vitro into adipogenic, chondrogenic, myogenic, and osteogenic cells in the presence of lineage-specific induction factors. In conclusion, the data support the hypothesis that a human lipoaspirate contains multipotent cells and may represent an alternative stem cell source to bone marrow-derived MSCs.

#### INTRODUCTION

THE THERAPEUTIC POTENTIAL of multilineage stem cells for applications such as tissue engineering and gene therapy is enormous. Conceptually, there are two general types of stem cells potentially useful for these applications: embryonic stem cells (ESCs) and autologous stem cells. Although theoretically appealing because of their pluripotentiality, the practical use of ESCs is limited due to potential problems of

<sup>1</sup>Dr. Zuk and Dr. Zhu are co-first authors.

<sup>2</sup>Laboratory for Regenerative Bioengineering and Repair, Departments of Surgery and Orthopaedic Surgery, UCLA School of Medicine, Los Angeles, California. <sup>3</sup>Division of Plastic and Reconstructive Surgery, University of Pittsburgh School of Medicine, Pittsburgh, Pennsyl-

vania.

#### What is Zuk Method?

#### Isolation and culture of stem cells-PLA and MSCs

Human adipose tissue was obtained from elective liposuction procedures under local anesthesia (HSPC #98-08 011-02). In this procedure, a hollow blunt-tipped cannula was introduced into the subcutaneous space through small (~1 cm) incisions. The cannula was attached to gentle suction and moved through the adipose compartment, mechanically disrupting the fat tissue. A solution of saline and the vasoconstrictor epinephrine was infused into the adipose compartment to minimize blood loss and contamination of the tissue by peripheral blood cells. The raw lipoaspirate (~300 cc) was processed according to established methodologies to obtain a stromal vascular fraction (SVF).<sup>33,34</sup> To isolate the SVF, lipoaspirates were washed extensively with equal volumes of phosphate-buffered saline (PBS), and the ECM was digested at 37°C for 30 min with 0.075% collagenase. Enzyme activity was neutralized with Dulbecco's modified Eagle's medium (DMEM), containing 10% FBS and centrifuged at  $1200 \times g$  for 10 min to obtain a high-density

SVF pellet. The pellet was resuspended in 160 mM NH<sub>4</sub>Cl and incubated at room temperature for 10 min to lyse contaminating red blood cells. The SVF was collected by centrifugation, as detailed above, filtered through a 100- $\mu$ m nylon mesh to remove cellular debris and incubated overnight at 37°C/5% CO<sub>2</sub> in control medium (DMEM, 10% FBS, 1% antibiotic/antimycotic solution). Following incubation, the plates were washed extensively with PBS to remove residual nonadherent red blood cells. The resulting cell population was termed a processed lipoaspirate (PLA), to distinguish it from the SVF obtained from excised adipose tissue. PLA cells were maintained at 37°C/5% CO<sub>2</sub> in noninductive control medium. Cells did not require specific FBS sera lots for expansion and differentiation (data not shown). For immunofluorescence studies, a population of MSCs was obtained from human bone marrow aspirates according to the protocol of Rickard *et al.*<sup>17</sup> and maintained in control medium. To prevent spontaneous differentiation, cells were maintained at subconfluent levels.

*Reference: Tissue Engineering. 7:212-228, 2001* 

- The Dr Zuk's adipose derived stem cell isolation technique was called Zuk Method in academics and a modified method from Zuk Method was also developed to satisfy human reinjection condition.
- For example, FBS (Fatal Bovine Serum) was replaced by Human Serum, and PBS (Phosphate Buffered Saline) was replaced by Normal Saline to satisfy human reinjection condition and this method was called Modified Zuk Method based on which some companies developed their ADSC isolation products (UNISTATION<sup>™</sup> also adopted this Modified Zuk Method).



#### What is **CAL** (Cell-Assisted Lipotransfer)?

- Since harvesting technique of ADSC was published by Dr Zuk (Tissue Engineering. 7:212-228, 2001), a number of research has been conducted for years to develop applications of ADSC.
- In 2008, Dr. Yoshimura, professor of plastic surgery department of Tokyo University in Japan, proved the efficacy of his theory that stem cells would increase the survival rate of fat transfer by angiogenesis / vascularization effect of stem cells, resulting in smoother oxygen supply and higher survival rate of injected fat (Aesth Plast Surg 2008;32:48-55).
- He conducted the research and he proved that the stem cell-rich lipotransfer can, in fact, enhance the fat survival rate up to 70%, compared to 30-40% of general lipotransfer.
- The research was a very extensive clinical test spending more than 2 years with breast augmentation cases of 40 patients (Aesth Plast Surg 2008;32:48-55) and with 6 facial treatment cases with LEP disease (Dermatol Surg 2008;34:1178–1185).
- And based on the success of the research, he named the stem cell rich fat transfer as CAL (Cell Assisted Lipotransfer).
- And many plastic surgeons, dermatologists, and urologists wanted to conduct the CAL surgery in their private practices because they had been tired of complaints from their customers for the high absorption rate of previous fat transfer.
- At that time, CAL was a very difficult surgery to be conducted by private practitioners though, due to the lack of facilities for stem cell isolation.
- However, after releases of several devices for ADSC isolation, including Multi Station, CHA-STATION, Celution, TGI, and etc, CAL treatment became available even for private practitioners. Thus CAL treatment became prevalent in private practices.

*Reference: Aesth Plast Surg 2008;32:48-55 Dermatol Surg 2008;34:1178–1185* 







#### Clinical Result of CAL by Dr. Yoshimura – Breast Augmentation



Before

After 24 months

Fig. 4 Clinical views of a patient in group A (patient 1). Preoperative (left) and postoperative (right) views at 24 months. A 22-year-old woman underwent breast augmentation using cell-assisted lipotransfer (CAL) (290 ml in each breast), with satisfactory results at 24 months. Her breast circumference increased by 5 cm, and her augmented breast mounds remained soft and natural appearing without injection scars or subcutaneous indurations

> Fig. 8 Clinical views of a patient in group C (patient 3). Preoperative (top) and postoperative (bottom) views at 24 months. A 30-year-old woman underwent breast augmentation with cell-assisted lipotransfer (CAL) (310 ml in each breast). Her breasts were dramatically augmented with an increase in breast circumference difference by 8 cm at 24 months. The breast mounds were soft with no subcutaneous

indurations. An original inframammary fold on the left breast is slightly visible, but injection scars are not visible



**Fig. 6** Clinical views of a patient in group C (patient 2). Preoperative (*top*) and postoperative (*bottom*) views at 12 months. A 32-year-old woman underwent breast augmentation with cell-assisted lipotransfer (CAL) (280 ml in each breast). Her breast circumference difference increased from 9 cm (baseline) to 14.5 cm (at 12 months). The breast mounds are soft and natural appearing with no visible injection scars

© 2008 by the American Society for Dermatologic Surgery, Inc. • Published by Wiley Periodicals, Inc. • ISSN: 1076-0512 • Dermatol Surg 2008;34:1178-1185 • DOI: 10.1111/j.1524-4725.2008.34256.x

Tokyo, Japan



	Clinical Result of CAL by Dr. Yoshimura							
Cell-Assisted Lipotransfer for F Clinical Use of Adipose-Derive	acial Lipoatrophy: Efficacy of d Stem Cells	TABLE 1. Summary of Pat	tients' Data					
Kotaro Yoshimura, MD,* Katsujiro Sato, MD,† Keita Inoue, MD,* Hirotaka Suga, MD,* Hiron Toshitsugu Hirohi, MD, <sup>§</sup> and Kiyonori Harii, N	Noriyuki Aoi, MD,* Masakazu Kurita, MD,‡ 11 Ero, MD,* Harunosuke Kato, MD,* 10‡		Method					
BACKGROUND Lipoinjection is a promising treatme remains to be established.	nt, but its efficacy in recontouring facial lipoatrophy	(	Non-CAL	)		CAL		
OBJECTIVE The objective was to evaluate the effice mentation of adipose-derived stem/stromal cells (A	acy and adverse effects of lipoinjection and supple- SCs) to adipose grafts.					$\smile$		
METHODS To overcome drawbacks of autologous called cell-assisted lipotransfer (CAL). In CAL, stro isolated from half of an aspirated fat sample and att the fat acting as a scaffold. This process converts rel performed conventional lipoinjection (non-CAL; <i>n</i> - ronby due to lunus perfundus or Parv-Romberg sy	lipoinjection, we have developed a novel strategy mal vascular fraction containing ASCs was freshly ached to the other half of aspirated fat sample with atively ASC-poor aspirated fat into ASC-rich fat. We 3) or CAL $(n=3)$ on six patients with facial lipoat- ndrome.	Case	1	2	3	1	2	3
RESULTS All patients obtained improvement in fac	ial contour, but the CAL group had a better clinical	Sex	Male	Female	Female	Male	Female	Female
significance $(p=.11)$ . Adipose necrosis was found corticosteroids.	in one non-CAL case who took perioperative oral	Age (vear)	25	42	42	35	33	48
CONCLUSION Our results suggest that CAL is both ventional lipoinjection for facial recontouring.	effective and safe and potentially superior to con-	Diagnosis		IED	IED	PRS	IED	IEP
The authors have indicated no significant interest with	commercial supporters.	Diagnosis		-		rno		LEF
		Severity grading	3	5	2	4	5	4
Facial lipoatrophy is a disfiguring and socially disabling problem that accompanies several in- herited and acquired diseases (see Ascher et al. <sup>1</sup> for review). Lupus erythematosus profundus and sele-	lipoatrophy but also bony defects. <sup>2</sup> Thus far there is no medical treatment to correct facial lipoatrophy.	Perioperative oral corticosteroids	-	+	-	-	+	-
phea) frequently give rise to facial lipoatrophy,	such as that associated with HIV infection with	Donor site	Abdomen	Abdomen	Abdomen	Abdomen	Abdomen	Abdomen
which is often the most problematic manifestation of the disease for the patients even though they require medication, such as oral corticosteroids, to suppress	moderate success, <sup>3,4</sup> while microsurgical tissue transfer has been a standard surgical treatment for recontouring facial defects seen in Parry-Romberg	Volume of injection (mL)	100	250	50	110	90	100
other symptoms associated with the systemic morbidities. Facial lipoatrophy is also seen in	syndrome. However, the surgical procedure leaves conspicuous scars on the face and the donor site. <sup>5,6</sup>	Follow-up periods (month)	10	9	9	13	9	9
patients with human immunodeficiency virus (HIV) infection and Parry-Romberg syndrome (hemifacial progressive atrophy or idiopathic	Autologous lipoinjection is a promising treatment for soft tissue augmentation because there is no associated incision scar or complications	Clinical results	Fair	Good	Fair	Excellent	Good	Good
hemifacial atrophy), which involves not only facial	associated with foreign materials. Although	Complications	_	Adipose necrosis	-	-	-	-
*Department of Plastic Surgery, University of Tokyo School Department of Plastic Surgery, Kyorin University School o Tokyo, Japan	of Medicine, Tokyo; <sup>†</sup> Cellport Clinic Yokohama, Yokohama; f Medicine, Tokyo; <sup>§</sup> Ritz Cosmetic Surgery Clinic Tokyo,							

LEP, lupus erythematosus profundus; PRS, Parry-Romberg syndrome.

**ADSC** enhances **Fat Survival Rate** up to **70%**, compared to 30% of previous Lipotransfer



















#### Limitations of Previous ADSC Devices

- Even though there have been several ADSC isolation devices, there have been also many technical problems, for example, air contamination possibility, long isolation time, high consumable costs, and so on.
- To solve those problems and to develop an upgraded form from the previous devices, NeoGenesis Co., Ltd. has conducted extensive research work for many years.

#### Development of **UNISTATION**<sup>TM</sup>

- NeoGenesis recruited research and marketing expertise in Adipose Derived Stem Cell field.
- NeoGenesis successfully developed a technique for preventing contamination possibility entirely.
- NeoGenesis successfully simplified the entire process of ADSC isolation with lowest consumable costs.
- NeoGenesis finally developed the simplest and most convenient ADSC isolation device ever, UNISTATION<sup>™</sup>.
- NeoGenesis added PRP, PRF isolation, and fat purification function to UNISTATION<sup>™</sup>.
- UNISTATION<sup>™</sup> is just now developed to satisfy the very needs of many surgeons who have yearned for treating patients with various cell therapies.

#### With UNISTATION<sup>TM</sup>, nothing is impossible!



## Detailed Information



#### Composition





#### Centrifuge & Incubator & Shaker in <u>One Device</u> <u>Smallest & Simplest system</u> for stem cell isolation

#### **Composition - Centrifuge**





16 UNIKITs can be installed. ( 4 UNIKITs per 1 Bucket)



Maximum capacity for centrifugation <u>50cc X 16 UNIKITs = 800cc High Volume Process at Once</u>

#### **Composition - Incubator**





#### Heat Release for Incubation from Silicon Heater



#### With Silicon Incubator installed inside, we keep the inner temp 37°C with the highest stability.

The most stable incubation system with **SILICON INCUBATOR** 

#### **Composition - Incubator**





#### The most stable incubation system with **SILICON INCUBATOR**

#### Composition - Shaker





The simplest Unified Shaking System

#### **Composition - Display**



#### Control Display

- $\bullet$  Easy Touch Screen  $\rightarrow$  Easy and Fast Order and Modification
- 10 Steps of Acceleration and Deceleration  $\rightarrow$  Suitable for Research with Sensitive Separated Layers
- RPM/RCF Auto-Calculation
- Timer Modifiable during Working
- Up to 1 Hour 59 Minutes 59 Seconds Time Setting
- Memory Function (100 Memories)
- Temperature Control (Room Temp ~ 40C)

#### Safety

- Auto Alarm System
- Imbalance Sensor
- Door Lock Safety
- Second Door Function
- RPM Error Sensor
- Motor Overheat Sensor

G	<b>()</b> h:m / m:s	J.C	1		Ш	SVF •	
			4	5	6	Fat Washing	
RPM TIME	TEMP ACC	PROG	7			PRP	<u> </u>
RCF	DEC	T HOU	+		-	PRF	((ا

#### The simplest **One Button Protocol** system & modern display



#### Waiting Display



#### Working Display: <u>Screen displays which process its working on</u>



- A: SVF in process
- B: Fat washing in process
- C: PRP in process
- D: PRF in process

\*\* Deceleration Step 0: Natural braking Steps 1~8: Slow stop Step 9: Fast stop



Specification	ns of UNISTATION™				
Max. RPF/RCF (Swing out)	3,600 rpm / 2,898 xg				
Max. Capacity (cc)	50cc x 4 x 4 (Total 800cc)				
Temp. range (℃)	Room Temperature ~ +40 ℃				
Pre-Warm Up Function	Yes (Set as A1 or A2)				
Time control	different for each mode				
RPM/RCF conversion	Yes				
Noise level (dB)	≤60				
Acc/Dec	9/10 steps				
Program memory	100				
Imbalance cutout	Yes				
Display	Blue LCD				
Safety lid lock	Yes				

#### Features

- Angle & Swing Rotor
- •Max 800ml (50ml UNIKITs x 4 holes x 4 buckets)
- •The Most Compact Design
- •Touch Screen Key Pad
- •RPM/RCF Auto-Calculation
- •100 Protocol Memorization
- •Separate Protocol Buttons for SVF, Fat, PRP, PRF each

#### **Consumables & Accessories**





UNIKIT (Kit body, cap, transfer.) \*\*Each part of the kit (UNIKIT, hand-piece, cap, transfer) are 1 set



We provide Consumable kits at lowest consumable price.

#### **Consumables & Accessories**





S750T-4B

TM96-4S

A100S-6

#### **Consumables & Accessories**



L												-											
S	500T	- 4B	No.12			B500	P	500 L i	d	TMQ6	-19	S	750T	- 4B	No.13			B750	Ē	3750 Li	d	TM96-	-4S
		Type Max. RPN Max. RCF Radius Dim.	Swin 4,500 4,392 194 r 262x	g Rotor ∠9 ) rpm 2xg mm 262x62	0°	Round Buck	et ity	Sealing Cap		ficroTiter P	late rack	Tube Capacity (ml) Dimensions Øx L mm	5 12x75	Type Max. RPM Max. RCF Radius Dim. 10 13x100	Swin 4,000 3,515 196.5 277x 15 16x114	g Rotor ∠90 0 rpm 5 x g 5 mm 277x55 50 29x104	P° F 7 250 61.8x123	Round Bucke (50ml capaci 500 69.5x170.2	et ty 750 97x152	Sealing Cap	15 17x120	licro Titer Pl	ate rack
Tube Capacity (ml)	5	10	15	50	100(85ml)	250	500	50	15	MTP		Tube Type			U	U	0	U		Ų	A INC.		
Dimensions Øx L mm	12x75	16x100	16x114	29x104	38x105.7	61.8x123	69.5x170.2	30x115	17x120				0	(Law)	den.	ad				at s	et a	-8	
Tube Type	П	П	Π	П	ñ	ñ			<b>6</b> ,0			Adapter		Û	Ű	Ũ	Ű				Ű	660	
	U		U	I U	U .	U		U		~		Tube per Adapter	26	26	19	7	1	1	-	5	14	5	-
				-			á					Tube per Rotor	104	104	76	28	4	4	4	20	56	20	4
			100	db				do	680			Max.RCF	3,433	3,417	3,438	3,288	3,463	3,458	3,515	3,452	3,413	3,335	2,925
Adapter		(22)	958 /	00	Imi	THE REAL	1	107	Mai			RPM	4,000	4,000	4,000	4,000	4,000	4,000	4,000	4,000	4,000	4,000	4,000
Adapter												boring Ø x L mm	13x70	14x89.5	17x99.5	29.2x81	62.2x99	75.5x99	99x105	29.5x99.5	17x89	27x65	-
	diffe	<b>u</b>	"	44	WIIIIP		$\sim$	-	All h			Color / Material	Magenta	Yellow	Green	White	Cyan	Violet	-	Blue	Sky Blue	Black	-
									-			Cat.NO.	T5-26	T10-26	T15-19	T50-7	T250-1	T500-1	-	T50-5C	T15-14C	TS50-5	-
Tube per Adapter	9	9	9	4	1	1	1	3	1	•		Max. length Ø x L	13x153.5	14x152.5	17x151	29.2x145	62.2x155	75.5x 55	99x158	29.5x146.5	17x152	27x90	-
Tube per Rotor	36	36	36	16	4	4	4	12	28	4			hand	mak		620							
Max.RCF	4,256	4,279	4,290	4,290	4,290	4,290	4,347	4,302	4,313	3,486		Adapter		F		P		P					
radius (mm)	188.0	189.0	189.5	189.5	189.5	189.5	192.0	190.0	190.5	154.0		Tube per Adapter	30	24	17	4	1	1					
RPM	4,500	4,500	4,500	4,500	4,500	4,500	4,500	4,500	4,500	4,000		Tube per Rotor	120	29	68	16	4	4					I
boring Ø x L mm	13x64	17x64	16.5x85	30x85.5	38.5x84.5	62x100	73x101	30x90	17x95.5			Max.RCF	3,336	3,336	3,336	3,336	3,336	3,336					
Color / Material	Croor	0.00	Khaki	Shy Blue	Brown	Viole*	Chartreuse	Rhua	Dink			radius (mm)	186.5	186.5	186.5	186.5	186.5	186.5					
Color / Material	Green	Cyan	NIIdKI	oky blue	brown	violet	Chartreuse	Diue	PILIK	•		RPM	4,000	4,000	4,000	4,000	4,000	4,000					i
Cat.NO.	TR5-9	TR10-9	TR15-9	TR50-4	TR100-1	TR250-1	TR500-1	TR50-3C	TR15-7C	•		boring Ø x L mm	13 x 55	13.5 x 80	17.3 x 80	29 x 80 White	62 x 80 Red	74 x 80 Blue					
Max. length Ø x L	13x150	17x151	16.5x151.5	30x151.5	38.5x151.5	62x151.5	73x154	30x152	17x152.5	-	1	Cat.NO.	T5-30	T10-24	T15-17	T50-4	T250-1	T500-1					
Combi 514R Max. RPM : 4,	500 <b>M</b> a	IX. RCF :	4,392 xg	9	C N	Combi 40 Max. RPN	0 <b>8</b> M∣: 4,000	Max. R	CF : 3,4	70 xg		Max. length Ø x L Combi 514R Max. RPM : 4,	13x148 000 <b>M</b> a	13.5x148	17.3x148 3,515 xg	29x148	62x148 C	74x148 Combi 40 Max. RPN	8 1:4,000	Max. R	CF : 3,5	15 xg	



## Protocols





#### Preparations for CAL surgery



UNISTATION<sup>™</sup> - 1 UNIKIT – 3 Closer – 1 40ml Suctioned fat – 1 20ml Collagenase – 1 5ml Patient's serum – 1 30ml Normal Saline - 1



#### <u>We provide the most</u> reasonable consumable price



#### Step 1 : Fat Washing (A1)

\*Before starting treatment, while suctioning fat, you may turn on the device and close the door and touch SVF button to set the device as A1 or A2. By this way, you can pre-warm up the device up to 37C for incubation process later. (Recommended)



1. Suctioned fat (40ml) and water for balance (40ml) in UNIKIT



2. Put them in UNISTATION<sup>™</sup> centrifuge



3. Touch SVF button once (Display shows A1)



4. Touch start button if A1 is shown



5. After centrifugation, remove RBC layer at the bottom and keep only pure fat



6. 20ml pure fat after removing RBC layer



#### Step 2 : Collagenase Digestion (A2)



1. Prepare 20ml of 0.1% Collagenase Solution



2. Transfer the Collagenase Solution into the pure fat

() htm/ma

2 3

6 0



3. Place Shaking Plate and put the pure fat with collagenase on it and close the door



6. Wait for 30 min for shaking incubation and take the fat out after finishing shaking incubation

# COMPARTING TO THE TANK OF THE

4. Touch SVF button once again (Display shows A2)

5. Touch Start button if A2 is shown



#### Step 3 : SVF Collection (A3)





4. After centrifugation, take out UNIKIT very slowly so that layers don't get mixed. (Stem cells are at the bottom within 5ml SVF)



5. Before removing cap, pull the hand piece upwards slightly so that stem cells doesn't stay in the cap. Slightly pull so that measurement marking stays under 1ml. \*The same should be done after centrifugation steps of A3 and A4.



6. Open the cap and use transfer to connect the Unikit to syringe to extract bottom 5ml of stem cell. Transfer the 5ml to new unikit.



#### Step 4 : Neutralization and washing of SVF (A4)







1. Transfer 5ml PPP into the 5ml SVF

2. Transfer 30ml normal saline into the 10ml of SVF and PPP

3. Put in UNISTATION™ centrifuge



4. Touch SVF button once again (Display shows A4)



5. Touch start button if A4 is shown



6. After centrifuation, transfer 5ml SVF at the bottom into 5ml or 10ml syringe



#### Preparation for injection



Transfer the final result of 5ml SVF into the 20ml pure fat for injection. Inject stem cell rich fat to the treatment area. (For IV injection, cell strainer must be used before injection.)



#### **Training videos**

#### Youtube Link of SVF / ADSC Isolation Video

http://youtu.be/Cpoz-8f\_qAw

Download Link of SVF / ADSC Isolation Video

https://docs.google.com/a/neogenesis.co.kr/file/d/0Bznemsb5skHkdkpESlVSeTMxNnc/edit

#### Youku Link of SVF / ADSC Isolation Video (Chinese)

http://v.youku.com/v show/id XODI0MjU2Njc2.html





ADSC for clinical use (SVF) -Fluorescent color

#### Purified ADSC after RBC lysis

#### SVF (Stem Cell) Pictures & Counting





**Fig. 9.** Multicolor flow cytometric analysis of cells in the bottom layer. CD45<sup>+</sup> cells were regarded as blood-derived cells. CD45<sup>-</sup> cells were regarded as adipose-derived cells and were processed to the next analysis. CD45<sup>-</sup>CD31<sup>-</sup>CD34<sup>+</sup> cells were regarded as adipose-derived stem/stromal cells (*ASCs*), whereas CD45<sup>-</sup>CD31<sup>+</sup>CD34<sup>+</sup> cells were regarded as endothelial cells.

\*Reference: Plast. Reconstr. Surg. 122: 111, 2008.

=60.7%

after RBC lysis

#### SVF (Stem Cell) Pictures & Counting



#### <u>Comparison of cell counting with other products</u>

20 00:57 AM	Details Sav	1999.4.23 08:32 Protocil: E6FALT
	Next sample	Save
	Total : 4.6x10 <sup>5</sup> mL	Screen Brightness Total : 3,5x10 /mL

Lot number	Yield (cells/ adipose mL)	Viability (%)	Adipose volume (mL)
1	241 800	85.0	50.0
2	449 500	90.0	50.0
3	278 471	88.0	121.0
4	189 000	80.6	55.0
5	380 586	89.7	145.0
6	231 701	86.0	97.0
Average	295 176	86.6	
SD	99 540	3.5	

#### <u>Comparison of cell counting with other products</u>

	UNISTATION <sup>TM</sup> $(n = 1)$	CHA-STATION (n = 2)	Multi Station $(n = 1)$	Cellution (n = 6)
SVF cells (per adipose 1cc)	<u>460,000</u>	300,000	350,000	295,176
ADSC cells (per adipose 1cc)	<u>287,274</u> (n=10)	180,000	210,000	177,105
Viability (%)	<u>90.1%</u>	85.1%	90.0%	86.6%

- 62.2% higher number of ADSC compared to Cytori
- 45.8% higher number of ADSC compared to CHA-Station
- 36.8% higher number of ADSC compared to Multi-Station



## Comparison & Reference



#### Comparison Study



	UNISTATION™ - NeoGenesis	CHA-STATION	Multi Station	Sceldis (HuriCell)	Celution-Cytori	Tissue Genesis Cell Isolation System-TGI
Separation Time	50 min	50 min	120 min	90 min	90 min	60 min
Automatic / Manual	Semi Automatic	Semi Automatic	Manual	Automatic	Automatic	Automatic
Closed Type (Contamination Prevention)	0	0	Х	0	0	0
Weight (Approximate)	70Kg	150Kg	200Kg	150Kg	70Kg	50Kg
Demo Possibility	0	Х	Х	Х	0	0
Cell Counting Function	Х	0	Х	Х	Х	Х
Max Volume	800cc	160cc	400cc	150cc	400cc	60cc
Price (Approximate)		USD 66,000	USD 40,000	USD 66,000	USD 150,000	USD 150,000
Consumable Price (One Set for One Time Operation, Approximate)	USD 150	USD 200	USD 100	USD 500	USD 2,000-11,000	USD 5,000
Good Things to Consider	Easy to Use / Low Consumable Price / Small Size / High Volume	Easy to Use / Cell Counting Function	High Volume / Low Consumable Price	Automatic	Automatic / High Volume	Automatic
Bad Things to Consider		Too Large / Small Volume / Low Durability	Contamination Possibility / Hard to Use / Too Large	Low Volume / Frequent A/S / Animal Enzyme Use / High Consumable Price	Expensive Consumable Price / Expensive Device Price / No Registration	Expensive Consumable Price / Expensive Device Price / No Registration / Low Volume
Customer Satisfaction	****	****	***	***	**	**

\*Red Colored means comparatively better than competitors



#### Semi-automated extraction and characterization of Stromal Vascular Fraction using a new medical device.

-Clinical Hemorheology and Microcirculation, vo. 64, no. 3, pp. 403-412, 2016

#### Abstract:

INTRODUCTION: A novel commercially available semi-automated device for the extraction of SVF promises sterility, consistent results and usability in the clinical routine. The aim of this work was to compare the quantity and quality of the SVF between the new system and an established manual laboratory method.

MATERIAL AND METHODS: SVF was extracted from lipo-aspirate both by a prototype of the semiautomated **<u>UNiStation</u><sup>™</sup>** (NeoGenesis, Seoul, Korea) and by hand preparation with common laboratory equipment. Cell composition of the SVF was characterized by multi-parametric flowcytometry (FACSCanto-II, BD Biosciences). The total cell number (quantity) of the SVF was determined as well the percentage of cells expressing the stem cell marker CD34, the leucocyte marker CD45 and the marker CD271 for highly proliferative stem cells (quality).

#### DISCUSSION: The semi-automated closed system provides a considerable amount of sterile SVF with high reproducibility.

Furthermore, the SVF extracted by both methods showed a similar cell composition which is in accordance with the data from literature. This semi-automated device offers an opportunity to take research and application of the SVF one step further to the clinic.

Semi-automated extraction and characterization of Stromal Vascular Fraction using a new medical device

Authors: Hanke, Alexander | Prantl, Lukas | Wenzel, Carina | Nerlich, Michael | Brockhoff, Gero | Loibl, Markus I Gehmert, Sebastian

#### Article Type: Research Article

Abstract: INTRODUCTION: The stem cell rich Stromal Vascular Fraction (SVF) can be harvested by processing lipo-aspirate or fat tissue with an enzymatic digestion followed by centrifugation. To date neither a standardised extraction method for SVF nor a generally admitted protocol for cell application in patients exists. A novel commercially available semi-automated device for the extraction of SVF promises sterility, consistent results and usability in the clinical routine. The aim of this work was to compare the quantity and quality of the SVF between the new system and an established manual laboratory method. MATERIAL AND METHODS: SVF was extracted from lipoaspirate both ... Show more

Keywords: Stromal vascular fraction, lipo-aspirate, adipose tissue, mesenchymal stem cells

DOI: 10.3233/CH-168124

Citation: Clinical Hemorheology and Microcirculation, vol. 64, no. 3, pp. 403-412, 2016

#### References



- (1) P. Zuk, et al. Multilineage Cells from Human Adipose Tissue: Implications for Cell-Based Therapies. TISSUE ENGINEERING 2001;7:212-228.
- (2) I. Hong, et al. Isolation and culture methods of adipose-derived stem cells. IFATS 2009 Oral Abstracts;17.
- (3) K. Yoshimura, et al. Cell-Assisted Lipotransfer for Cosmetic Breast Augmentation: Supportive Use of Adipose-Derived Stem/Stromal Cells. Aesth Plast Surg 2008;32:48-55.
- (4) K. Yoshimura, et al. Cell-Assisted Lipotransfer for Facial Lipoatrophy: Efficacy of Clinical Use of Adipose-Derived Stem Cells. Dermatol Surg 2008;34:1178–1185.
- (5) H. Suga, et al. Numerical Measurement of Viable and Nonviable Adipocytes and Other Cellular Components in Aspirated Fat Tissue. *Plast. Reconstr. Surg.* 2008;122:113.
- (6) J. Mitchell, et al. Immunophenotype of Human Adipose-Derived Cells: Temporal Changes in Stromal-Associated and Stem Cell–Associated Markers. STEM CELLS *2006;24:380.*
- (7) A Hanke et al. Semi-automated extraction and characterization of Stromal Vascular Fraction using a new medical device. Clin Hemorheol Microcirc 64 (3), 403-412. 2016.



## NEOGENESIS •

#### Global Leader Enterprise of Autologous Cell Therapies

Address. #302, 146, Seonyu-ro, E&C Dream Tower, Yeongdeungpo-gu, Seoul, South Korea (07255) Tel. +82-2-583-1348 Fax. +82-2-2038-0848 Email. info@neogenesis.co.kr Homepage. www.neogenesis.co.kr



## LUNASTEM<sup>TM</sup> Dual Fluorescence Cell Counter





æ

Advanced Dual fluorescence optics

#### <u>Fast</u>

Images generated from three channels in 30-60 s

Accurate Distinguishes nucleated and non-nucleated cells

#### **Convenient**

East-to-use software with a simple counting procedure

#### **LUNA**STEM<sup>™</sup>

The LUNA-STEM<sup>™</sup> is a stand-alone compact instrument that combines dual fluorescence optics and image analysis software. Its interactive touchscreen makes accurate cell counting a simple experience.







#### **Dual Fluorescence + Bright Field Optics**

The LUNA-STEM<sup>™</sup> is a quantum leap forward for automated cell counting and viability analysis. Using Acridine Orange (AO) and Propidium Iodide (PI) to stain live and dead cells, respectively, the LUNA-STEM<sup>™</sup> provides highly sensitive and accurate results for most cell types, including stem cells and SVF cells.

The various cells within the SVF can be recognized and counted as non-nucleated cells or live and dead nucleated cells, which is not possible with other cell counters.





LUNASTEM"	E

#### **PhotonSlide**<sup>™</sup>

The PhotonSlide<sup>™</sup> has been developed to have ultra-low autofluorescence. Its patented precision chamber has a height that allows for the even distribution of cells throughout the counting chamber.







#### **Nucleated & Non-Nucleated Cell Analysis**





LUNA-STEM <sup>®</sup>	







#### **Interactive Softwar Interface**

#### **Powerful On-Board Analysis**

Built-in software automatically generates cell viability data upon cell counting. For validation purposes, live and dead cells are tagged with green and red circles, respectively.

#### **Image Overlay**

Images are captured from three channels (brightfield, green, and red) and can be merged directly on the screen. The brightness of each color can be adjusted individually for accurate detection. All images can be saved to an external USB drive.

#### **Data Report**

Data can be exported as a PDF report or CSV file. The LUNA-STEM<sup>™</sup> generates comprehensive PDF reports complete with the cell count and viability data, cell images, and relevant histograms.

#### Luna<sup>™</sup> Printer

The results of each count may be printed out immediately with the LUNA<sup>™</sup> printer.





#### **Specifications**

Sample Volume	10 µL
Cell Counting Time	30-60 sec (depending on sample conditions)
Cell Concentration Range	5×10 <sup>4</sup> - 2×10 <sup>7</sup> cells/mL(optimal range)
Cell Size Range	Detectable Range: 1 - 90 μm Optimal Range: 5 - 60 μm
Excitation wavelength	470 ±20 nm
Emission wavelength	530 ±25 nm, 600 nm (LP)
Light Source	LED
Image Resolution	5 MP
LCD Display	7 inches (800 x 480 pixels)
Dimensions (W×D×H)	22 × 21 × 9 cm (8.6 × 8.3 × 3.5 inch)
Weight	1.8 kg (4 lb) *without the AC adaptor

#### **Ordering Information**

Cat #	Product	Quantity
L30001	LUNA-STEM™ Automated Fluorescence Cell Counter	1 unit
L12005	PhotonSlide™, 50 Slides	1 box
L12006	PhotonSlide™, 500 Slides	10 box
L12007	PhotonSlide™, 1000 Slides	20 box
F23102	LUNA™ Fluorescence Calibration Beads	1 x 0.5 mL
F23001	Acridine Orange/Propidium lodide (AO/PI) Stain	2 x 0.5 mL
F23002	Acridine Orange (AO) Stain	2 x 0.5 mL
F23003	Propidium lodide (PI) Stain	2 x 0.5 mL
F23212	Cell Dilution Buffer	5 x 20 mL
P10001	LUNA™ Printer	1 unit
P10002	LUNA™ Printer Paper - thermal, 700 prints	3 x 2 rolls
P13001	LUNA™ Printer Cleaning Pen	1 unit
U10004	USB Drive, 4 GB	1 unit



## NEOGENESIS

Global Leader Enterprise of Autologous Cell Therapies

Address. #302, 146, Seonyu-ro, E&C Dream Tower, Yeongdeungpo-gu, Seoul, South Korea (07255) Tel. +82-2-583-1348 Fax. +82-2-2038-0848 Email. info@neogenesis.co.kr Homepage. www.neogenesis.co.kr

## **Genesis PRP**

The Simplest and Easiest Way to Collect the Most Concentrated PRP (Platelet Rich Plasma)



NeoGenesis Co., Ltd. Copyright 💿 All Rights Reserved

#### **Special Features of Genesis PRP**

- Transparent cylindrical kit allows a clear visibility and easier extraction
- The **curved neck** design reduces the possible cell loss
- PRP and PPP can be extracted from a single kit by one-time centrifugation
- High concentration rate of Platelets (about 13 times)
- Does not use a needle when collecting PRP or PPP
- Can control the volume of Plasma according to the treatment plan
- Closed system to prevent from air contamination possibility
- Always yields the maximum amount or PRP regardless of the operator's skill
- Compatible with any swing type centrifuge
- Competitive price



Genesis PRP

#### **Instruction Manual of Genesis PRP – Preparation**



The following accessories do not come with Genesis PRP kits.
The size of syringes can vary according to the user's preference.

- Needle (3ea)
- Anti-coagulant (1ea)
- 3cc syringe for PRP extraction (1ea)
- 10cc syringe for PPP extraction (1ea)
- 20cc syringe for drawing blood (1ea)
- Genesis PRP kit (1ea)
- A tube for balancing Centrifuge (1ea)
- Buffy controller (1ea)



1) Collect **1.5cc of anti-coagulant** into a syringe



2) Using a different needle, draw **15cc of patient's blood** into the same syringe and gently invert the syringe.

#### Genesis PRP

#### Instruction Manual of Genesis PRP – PRP & PPP Extraction



1) Inject the prepared blood into the kit by using an 18 gauge needle.



2) Place another Genesis PRP kit (or any kind of tube is fine) filled with water/saline in the same weight to balance the centrifuge .



3) Centrifuge the kits at 1700G (RCF) for 5 minutes, or press the "PRP" button in UNiStation; it will automatically centrifuge the kits at 3,000 rpm for 5 minutes.



4) Results divided into 3 different layers



5) Replace the bottom cap of the kit with Buffy controller's cap, and assemble the pusher.



6) Turn the bottom of the Buffy controller counterclockwise until Plasma reaches the top.

#### Genesis PRP

#### Instruction Manual of Genesis PRP – PRP & PPP Extraction



7) Connect a PPP collecting syringe, and rotate the controller until the RBC layer reaches the "1.0" line on the kit.



8) Connect a **PRP** collecting syringe, and rotate the pusher until the RBC layer reaches the top.



9) PPP (left) and PRP (right) extraction finished.



10) PRP activation is possible with NeoGenesis' PRP Activation Kit



If PPP extraction is not necessary, keep turning the controller without connecting a syringe so that Plasma can overflow into the empty space, and then connect the PRP collecting syringe and extract as much PRP as needed.



#### Youtube Link of Genesis PRP Manual Video

http://youtu.be/494YeLv6mZE

#### Download Link of Genesis PRP Manual Video

https://drive.google.com/file/d/oBwNSediveEuBbop6MWhoNGpxVnM/edit?usp=sharing

NeoGenesis Co., Ltd. Copyright 💿 All Rights Reserved.



Global Leader Enterprise of Autologous Cell Therapies

Address: #302, 146, Seonyu-ro, E&C Dream Tower, Yeongdeungpo-gu, Seoul, South Korea (07255) Tel: +82-2-583-1348 Fax: +82-2-2038-0848 Email: info@neogenesis.co.kr Homepage: www.neogenesis.co.kr