



Miltenyi Biotec

# MACS<sup>®</sup> Cell Separation

Select the best



# MACS® Technology

## One portfolio for all your cell isolation needs

MACS® Technology enables the magnetic separation of cell populations based on surface antigens. It is a fast and gentle method for the isolation of viable and functional cells by labeling epitopes with specific antibodies conjugated to superparamagnetic beads.

The MACS Technology portfolio provides a broad range of options for the isolation of virtually any cell type. Thereby, you enjoy the freedom to choose the cell isolation method that is best for your cells and your specific requirements. Our portfolio offers consistent, reliable cell separation solutions across basic and clinical research.

With MACS Technology, you are sure to select the best.

### Starting material

		Single-cell suspensions e.g. PBMCs, dissociated tissues (incl. tumors)	Blood products e.g. whole blood, buffy coat, apheresis products
Isolation strategy	Positive selection	<b>MACS® MicroBeads</b> <ul style="list-style-type: none"><li>• Columns</li><li>• Nano-sized MicroBeads</li><li>• UltraPure and REAlease™ MicroBeads</li></ul>	<b>StraightFrom® MicroBeads</b> <ul style="list-style-type: none"><li>• Columns</li><li>• Nano-sized MicroBeads</li></ul>
	Untouched isolation	<b>MACS Cell Isolation Kits</b> <ul style="list-style-type: none"><li>• Columns</li><li>• Nano-sized MicroBeads</li></ul>	<b>MACSxpress® Beads</b> <ul style="list-style-type: none"><li>• Column-free</li><li>• Micro-sized MACSxpress Beads</li></ul>

## MACS® MicroBeads

See what makes MACS Technology the most-cited cell isolation technology.

[See page 4](#)



## MACS® Columns

Learn about the advantages of MACS Columns – tested and trusted.

[See page 5](#)



## The MACS® Technology advantage

Discover the advantage of using MACS MicroBeads and Columns.

[See page 6](#)



## Cell separation

- Cell isolation from single-cell suspensions and dissociated tissues
- Cell isolation directly from blood products
- The next step in flexibility – label-free cells and challenging samples

[See pages 7–9](#)



## Manual and automated cell isolation

Choose the right cell isolation method for your specific needs.

[See page 10](#)




# MACS® MicroBeads

## MACS® MicroBeads – proven technology for basic research and clinical applications

MACS® MicroBeads are 50-nm superparamagnetic particles that are conjugated to highly specific antibodies against a particular cell surface antigen. Due to their small size, the beads do not activate cells. Furthermore, MACS MicroBeads do not have to be removed for any downstream application.


- MACS MicroBead Technology gives you the most flexible, most proven method for cell separation
- Minimal cell labeling with nano-sized MicroBeads ensures preservation of cellular integrity and characteristics
- Used in over 55,000 clinical cellular treatments to date



**LEARN MORE** 

For more information on MACS MicroBeads please visit  
▶ [miltenyibiotec.com/microbeads](https://miltenyibiotec.com/microbeads)

MACS MicroBead Technology owes its longstanding success to the ingenious combination of nano-sized superparamagnetic beads and a strong magnetic field in our MACS Columns. Only this technology ensures minimal labeling of target cells and the preservation of cellular properties. Cell separation with MACS MicroBeads is based on three easy steps: magnetic labeling, magnetic separation, and elution of labeled cells (fig. 1).



**Figure 1:** It only takes three easy steps to get viable cells with high yield and purity from your sample.

**VIDEO** 

Watch how to isolate cells in three easy steps at  
▶ [miltenyibiotec.com/3-easy-steps](https://miltenyibiotec.com/3-easy-steps)

# MACS® Columns

## MACS Columns – maximal magnetic power for minimal cell labeling

At the heart of MACS® MicroBead Technology is the MACS Column, containing a matrix composed of ferromagnetic spheres covered with a cell-friendly coating.



**Figure 2:** MACS Columns were developed for the fast separation of any cell type labeled with MACS MicroBeads.

When the column is placed in a MACS Separator, the spheres amplify the magnetic field by 10,000-fold, thus inducing a strong magnetic force within the column. The magnetic field efficiently retains cells labeled with the small, nano-sized beads.




**Figure 3:** MACS Column placed in a MidiMACS™ Separator.

[LEARN MORE](#)


Tailored formats for excellent results – find the optimal column for your cells at [▶ miltenyibiotec.com/columns](https://www.miltenyibiotec.com/columns)

The spacious matrix inside the MACS Columns ensures that unlabeled cells can easily flow through while minimally labeled cells (fig. 4) are gently yet effectively retained (fig. 5). This minimizes stress on the cells and allows for efficient washing while preventing cell aggregation.



**Figure 4:** Without the use of a MACS Column, extensive labeling or large beads are needed for an adequate magnetic retention. Only when using MACS Columns, the amplification of the magnetic force ensures effective cell retention with minimal labeling using the small beads.

## MACS® Columns enable gentle flow of cells. No pressure, sticking, or compression.



### Column matrix


- Designed and manufactured to deliver maximum purity and yield of viable cells
- The space between the spheres is about 20 times the size of lymphocytes
- Gentle to cells as they can freely flow through

**Figure 5:** The MACS Column at a glance. Cells move freely between the spheres inside the column and are only retained by magnetic forces.

# The MACS® Technology advantage

## Select the best by combining MACS® Columns and MicroBeads


### Advantage of column-based technology



**MACS® MicroBead Technology**


Strong magnetic force

---



**Minimal labeling suffices**


---



**Benefits of minimal labeling**

- No non-specific labeling
- No cell activation
- No alteration of cell characteristics


### Disadvantages of column-free technology



**Competing column-free technology**


Weak magnetic force

---



**Massive labeling required**

---



**Consequences of massive labeling**

- Non-specific labeling
- Cell activation
- Alteration of cell characteristics


**Figure 6:** Human PBMCs were either labeled with MACS CD3 MicroBeads for the isolation of T cells with a MACS Column or with other nano-sized beads for column-free isolation of the same cell type. Scanning electron microscopy showed (A) no visible labeling on the cell surface after isolation with MACS MicroBeads and MACS Columns, whereas (B) excessive labeling became obvious (indicated by arrows) after isolation with column-free technology from another manufacturer.

### YOUR BENEFITS



Why you select the best with MACS Technology:

- Effective separation: maximum purity and recovery
- Small bead size and minimal labeling: preserved cell functionality
- No cell stress: highest cell viability
- Free epitopes, no bead aggregation, no epitope cross-linking: full downstream compatibility



**Figure 7:** Light microscopic analysis of human PBMC cultures labeled with MACS CD3 MicroBeads or with nano-sized beads from another manufacturer. (A) No bead accumulation in cell culture observed with MACS MicroBeads. (B) Clearly visible bead aggregation (brown) with the other technology.


# Cell isolation from single-cell suspensions and dissociated tissues

## MACS® MicroBeads and MicroBead Kits

### Straightforward positive selection of target cells based on specific markers

The strong magnetic field generated by the matrix in the MACS® Column allows for minimal labeling of target cells with nano-sized MicroBeads. This ensures that plenty of surface epitopes remain free for subsequent fluorescent staining and flow cytometry analysis. Moreover, low labeling concentrations and the small size of MACS MicroBeads do not lead to activation of target cells (fig. 8).

- The least manipulative positive selection method
- Preservation of cell functionality due to optimal labeling
- Biodegradable: labeled cells are ready for downstream applications




**Figure 8:** Human B cells were enriched using MACS CD19 MicroBeads or a column-free positive selection method from another manufacturer. Subsequently, cells were cultured for 7 days in the presence or absence of the B cell stimulation reagents CD40-Ligand/Anti-His antibody and IL-4. Activation markers (CD69, CD80, and CD86) were measured by flow cytometry directly after cell isolation and after cultivation with and without stimulation. MACS MicroBeads did not alter the status of the target cells, whereas the column-free method led to the activation of B cells in the absence of stimulation reagents.

## MACS® Cell Isolation Kits


### Depletion of non-target cells to obtain pure, truly untouched cells

MACS® Cell Isolation Kits contain a cocktail of titrated antibodies and MACS MicroBeads for indirect magnetic labeling (fig. 9). They are the preferred choice if binding of antibodies to the target cells is not desired. Minimal labeling of unwanted cells with MACS MicroBeads avoids non-specific labeling of target cells, leaving the target cells truly untouched (fig. 10). In contrast, column-free methods based on nano-sized beads from other manufacturers require high concentrations of labeling reagents resulting in non-specific labeling of the target cell fraction.

- High purity and recovery rates
- Fully compatible with any downstream application
- No non-specific labeling of target cells



**Figure 9:** Non-target cells are magnetically labeled and depleted. During separation, the unlabeled target cell type is collected in the flow-through fraction. The labeled non-target cells are retained within the column.



**Figure 10:** Monocytes were enriched by depletion of unwanted cells using (A) the MACS Monocyte Isolation Kit II, human or (B) a column-free kit for human monocyte isolation from another manufacturer. Staining of monocytes (red) and nano-sized beads (green) after isolation showed non-specific labeling of the target cells when using column-free kits, while MACS Technology provided truly untouched cells.

# Cell isolation directly from blood products

## StraightFrom® Technology

### Cell isolation directly from blood products without density gradient centrifugation

StraightFrom® MicroBeads allow magnetic isolation of various leukocyte subsets from different starting materials by positive selection. With these kits, isolation of leukocyte subsets has never been easier and quicker. In contrast to conventional methods, StraightFrom Technology does not require density gradient centrifugation (fig. 11).

- Start directly with whole blood, buffy coat, and leukocyte reduction system chamber (LRSC)
- The isolated target cells are immediately ready for any downstream application
- Simple protocol with only a few handling steps

StraightFrom MicroBeads protocol	Conventional protocol
Transfer and dilute sample	Transfer and dilute sample
Magnetic labeling	Layer density gradient medium
Cell separation	Density gradient centrifugation without brake
<b>&lt;30 min</b>	PBMC separation and wash
	Cell count
	Magnetic labeling
	Cell separation
	<b>&gt;2 h</b>


**Figure 11:** Comparison of the StraightFrom MicroBeads protocol with conventional protocols, demonstrating the simplicity and short hands-on time.

## MACSxpress® Technology

### With high speed to untouched target cells

MACSxpress® Technology enables the fastest large-scale isolation of untouched cells directly from whole blood – without the need for any centrifugation. Micro-sized MACSxpress Beads allow for minimal labeling to prevent non-specific labeling and activation of target cells. Non-target cells are removed by immunomagnetic depletion. Simultaneously, erythrocytes are sedimented to yield target cells of exceptional purity (fig. 12).

- Go from whole blood to pure cells within 20 minutes
- Obtain untouched target cells directly from whole blood
- No density gradient centrifugation, erythrocyte lysis or cell counting required



**Figure 12:** MACSxpress Technology allows the isolation of cells from whole blood within 20 minutes.




# The next step in flexibility – label-free cells and challenging samples

## REAl ease™ Technology

### Get bead- and label-free cells

REAl ease™ MicroBead Kits have been developed for positive selection of target cells from PBMCs. REAl ease MicroBead Technology relies on recombinantly engineered antibody fragments instead of antibodies to label specific cell surface markers. The antibody fragments have a low affinity for cell surface epitopes. However, when the fragments are multimerized as a complex, they bind epitopes with high avidity and enable effective magnetic cell separation. REAl ease Technology controls the multimer/monomer state of the fragments and thus triggers the release of monomerized antibody fragments from the cell surface after isolation. Ultimately, the isolated cells are free from antibody fragments and magnetic labels.

- Bead-free cells: suited for second round of magnetic labeling
- Label-free cells: the epitope of a marker becomes completely available again
- Recombinantly produced: lot-to-lot consistency allows for reproducible results



LEARN MORE 


Learn more about REAl ease MicroBead Technology at [▶ miltenyibiotec.com/release-microbeads](https://miltenyibiotec.com/release-microbeads)

## UltraPure MicroBeads

### Minimize debris for high-quality results

UltraPure MicroBeads have been particularly optimized for use with challenging samples. The unique formulation provides compelling benefits particularly when starting with materials that contain large amounts of cell debris or low numbers of target cells. UltraPure MicroBeads greatly improve recovery and purity of the sorted population by specifically enriching viable target cells (fig. 13).

- Optimized formulation to minimize debris
- High cell purity, even from challenging starting materials
- As easy to use as the classic MACS® MicroBeads



**Figure 13:** CD34<sup>+</sup> cells were isolated with the column-based CD34 MicroBead Kit UltraPure (upper plots) or with a column-free positive selection method from another manufacturer (lower plots). The cell population purified with MACS MicroBeads UltraPure showed greatly reduced amounts of debris compared to the column-free method.

LEARN MORE 

Find the separation strategy that best fits your needs at [▶ miltenyibiotec.com/separation-strategies](https://miltenyibiotec.com/separation-strategies)

# From manual to fully automated high-throughput cell isolation

## Manual separation

Ease-of-use with manual MACS® Separators for simple and straightforward setups in any lab.

- The ideal solution for low-throughput experiments
- Proven technology in over 30,000 publications
- Perfectly tailored solutions for your experimental needs



Figure 14: Manual MACS Separators equipped with MACS Columns.

LEARN MORE 

First steps into MACS Technology – manual MACS Separators at a glance  
► [miltenyibiotec.com/separators](https://miltenyibiotec.com/separators)

## autoMACS® Pro Separator

Fully automated benchtop instrument for magnetic cell separation of multiple samples.

- Walk-away automation with cell labeling and isolation of up to six samples
- Standardized cell separation for reproducible, user-independent results
- Intuitive, easy-to-use software interface for a multi-user environment



Figure 15: Fully automated labeling and separation for the most convenient way to obtain pure cell populations with the autoMACS Pro Separator.

LEARN MORE 

Explore the autoMACS® Pro Separator features and watch the video at  
► [miltenyibiotec.com/automacs](https://miltenyibiotec.com/automacs)



SUPPORT 

Miltenyi Biotec offers comprehensive technical support for both new and advanced users alike. Our experienced technical support teams have the knowledge and expertise to answer your questions.

You can reach us at your convenience by e-mail, phone, or online in our forums and Live Chat – find all the information at  
► [miltenyibiotec.com/support](https://miltenyibiotec.com/support)

No question is too big or small.


## MultiMACS™ Cell24 Separator Plus

Efficient, semi-automatic cell isolation of large sample volumes or numbers.

- Convenient and easy handling of up to 24 samples in parallel or large sample volumes
- Compatible with any starting material and cell separation strategy
- Reliable, standardized process for reproducible results



**Figure 16:** Functional design for the isolation of large sample numbers or volumes with the semi-automated MultiMACS Cell24 Separator Plus.

**LEARN MORE** 

Simultaneous multisample magnetic cell separation with the MultiMACS™ Cell24 Separator Plus  
▶ [miltenyibiotec.com/multimacs](https://miltenyibiotec.com/multimacs)

## MultiMACS™ X

Walk-away solution for high-throughput setups – the next level of automated cell separation.

- The benefits of the MultiMACS™ Cell24 Separator Plus integrated into a liquid handler for minimal hands-on time
- Tailored solutions for your specific application
- Sample tracking, run reports, and LIMS integration



**Figure 17:** Full automation, high-throughput processing, and sample tracking for true walk-away cell isolation with the MultiMACS X.

**LEARN MORE** 

MultiMACS X – designed to speed up automated cell separation  
▶ [miltenyibiotec.com/multimacsx](https://miltenyibiotec.com/multimacsx)

► [miltenyibiotec.com/cellseparation](http://miltenyibiotec.com/cellseparation)



**Miltenyi Biotec**

**Germany/Austria/  
Switzerland**

Miltenyi Biotec GmbH  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany  
Phone +49 2204 8306-0  
Fax +49 2204 85197  
[macs@miltenyibiotec.de](mailto:macs@miltenyibiotec.de)

**USA/Canada**

Miltenyi Biotec Inc.  
2303 Lindbergh Street  
Auburn, CA 95602, USA  
Phone 800 FOR MACS  
Phone +1 530 888 8871  
Fax +1 877 591 1060  
[macs@miltenyibiotec.com](mailto:macs@miltenyibiotec.com)

**Australia**

Miltenyi Biotec  
Australia Pty. Ltd.  
Unit 16A, 2 Eden Park Drive  
Macquarie Park NSW 2113  
Australia  
Phone +61 2 8877 7400  
Fax +61 2 9889 5044  
[macs@miltenyibiotec.com.au](mailto:macs@miltenyibiotec.com.au)

**Benelux**

Miltenyi Biotec B.V.  
Schipholweg 68 H  
2316 XE Leiden  
The Netherlands  
[macs@miltenyibiotec.nl](mailto:macs@miltenyibiotec.nl)  
**Customer service**

**The Netherlands**  
Phone 0800 4020120  
Fax 0800 4020100

**Customer service Belgium**

Phone 800 94016  
Fax 0800 99626

**Customer service Luxembourg**

Phone 800 24971  
Fax 800 24984

**China**

Miltenyi Biotec Technology &  
Trading (Shanghai) Co., Ltd.  
Rooms 2303 and 2309  
No. 319, Xianxia Road  
Changning District  
200051 Shanghai, P.R. China  
Phone +86 21 62351005  
Fax +86 21 62350953  
[macs@miltenyibiotec.com.cn](mailto:macs@miltenyibiotec.com.cn)

**France**

Miltenyi Biotec SAS  
10 rue Mercoeur  
75011 Paris, France  
Phone +33 1 56 98 16 16  
Fax +33 1 56 98 16 17  
[macs@miltenyibiotec.fr](mailto:macs@miltenyibiotec.fr)

**Italy**

Miltenyi Biotec S.r.l.  
Via Paolo Nanni Costa, 30  
40133 Bologna  
Italy  
Phone +39 051 6 460 411  
Fax +39 051 6 460 499  
[macs@miltenyibiotec.it](mailto:macs@miltenyibiotec.it)

**Japan**

Miltenyi Biotec K.K.  
Nittsu-Eitai Building 5F  
16-10 Fuyuki, Koto-ku,  
Tokyo 135-0041, Japan  
Phone +81 3 5646 8910  
Fax +81 3 5646 8911  
[macs@miltenyibiotec.jp](mailto:macs@miltenyibiotec.jp)

**Nordics and Baltics**

Miltenyi Biotec Norden AB  
Scheelevägen 17  
223 70 Lund  
Sweden  
[macs@miltenyibiotec.se](mailto:macs@miltenyibiotec.se)

**Customer service Sweden**

Phone 0200-111 800  
Fax 046-280 72 99

**Customer service Denmark**

Phone 80 20 30 10  
Fax +46 46 280 72 99

**Customer service**

**Norway, Finland, Iceland,  
and Baltic countries**  
Phone +46 46 280 72 80  
Fax +46 46 280 72 99

**Singapore**

Miltenyi Biotec Asia Pacific Pte Ltd.  
100 Beach Road  
#28-06 to 28-08 Shaw Tower  
Singapore 189702  
Phone +65 6238 8183  
Fax +65 6238 0302  
[macs@miltenyibiotec.com.sg](mailto:macs@miltenyibiotec.com.sg)

**South Korea**

Miltenyi Biotec Korea Co., Ltd  
Arigi Bldg. 8F  
562 Nonhyeon-ro  
Gangnam-gu  
Seoul 06136, South Korea  
Phone +82 2 555 1988  
Fax +82 2 555 8890  
[macs@miltenyibiotec.co.kr](mailto:macs@miltenyibiotec.co.kr)

**Spain**

Miltenyi Biotec S.L.  
C/Luis Buñuel 2  
Ciudad de la Imagen  
28223 Pozuelo de Alarcón (Madrid)  
Spain  
Phone +34 91 512 12 90  
Fax +34 91 512 12 91  
[macs@miltenyibiotec.es](mailto:macs@miltenyibiotec.es)

**United Kingdom**

Miltenyi Biotec Ltd.  
Almac House, Church Lane  
Bisley, Surrey GU24 9DR, UK  
Phone +44 1483 799 800  
Fax +44 1483 799 811  
[macs@miltenyibiotec.co.uk](mailto:macs@miltenyibiotec.co.uk)

[www.miltenyibiotec.com](http://www.miltenyibiotec.com)

Miltenyi Biotec provides products and services worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. autoMACS, MACS, the MACS logo, MACSxpress, MidiMACS, MultiMACS, REAlease, and StraightFrom are registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. Copyright © 2017 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.